

GPC reader is an R Shiny application for routine calculations of M_p , average molecular weights (M_n and M_w), and PDI from GPC chromatograms (Gel permeation chromatography/size exclusion chromatography (SEC)). You can upload multiple files, inducing calibration files, and save results as csv file.

All calculations are based on the formulas given in [1].

$$M_n = \frac{\sum H_i}{\sum \frac{H_i}{M_i}} \quad M_w = \frac{\sum H_i M_i}{\sum H_i} \quad PDI = \frac{M_w}{M_n}$$

M_n number-average molecular weight;

M_w weight-average molecular weight;

PDI polydispersity index;

H_i height at the i-th point (detector signal intensity);

M_i molecular weight at the i-th point.

M_p is calculated as

$$M_p = 10^{(t * slope + intercept)} \quad \text{for a linear calibration curve}$$

$$M_p = 10^{(t^3 * coef\ 3 + t^2 * coef\ 2 + t * coef\ 1 + intercept)} \quad \text{for a polynomial calibration curve}$$

t is retention time (or retention volume).

App plots the molecular weight distribution (MWD), specifically differential MWD.

$dW/d(\log M_i)$ is calculated as

$$\frac{dW_i}{d(\log M_i)} = \frac{dW_i}{dV_i} \frac{dV_i}{d(\log M_i)}$$

Where

$$\frac{dW_i}{dV_i} = \frac{H_i}{\sum \frac{H_i}{\Delta V}}$$

And ΔV is constant time (volume) interval.

Sidebar and Tabs

The application has a sidebar and two tabs.

The sidebar panel and the first tab "Calibration" are intended for setting the coefficients of the calibration curve. On the first tab, you will have two tables and two plots:

- 1) **Table 1. Calibration samples.** This table shows the names of the corresponding files, Mp values and peak retention times for each calibration standard.
- 2) **Plot 1. GPC chromatogram of standards.** Detector signal intensity versus time.
- 3) **Table 2. Calibration coefficients.** It shows the calibration plot coefficients calculated from calibration files or entered in the sidebar.
- 4) **Plot 2. GPC calibration curve.** The logarithm of the molecular weight of polymer standards, vs the retention time (the elution volume).

GPC reader

Calibration:
☒ linear
☐ polynomial 3

Enter the Intercept

and the Slope

for a polynomial calibration curve, enter the 2nd

and 3rd coefficients

☐ use these coefficients

☐ use these coefficients

System peak, calibration

System peak, samples

☐ Add tails

Change the peak boundaries:

☐ shift start time

☐ shift end time

CalibrationData

Select calibration file(s)
Browse... 7 files

Table 1. Calibration samples

	FileName	Mp	RT
1	PMMA200kDa.txt	200000	8.21666666666667
2	PMMA102kDa.txt	102000	8.48333333333333
3	PMMA31_8kDa.txt	31800	9.18333333333333
4	PMMA23200Da.txt	23200	9.45
5	PMMA12_8kDa.txt	12800	10.35
6	PMMA5_93kDa.txt	5930	11.0166666666667
7	PMMA2_13kDa.txt	2130	11.45

Showing 1 to 7 of 7 entries

Previous1Next

☒ one sample per file

Showing 1 to 7 of 7 entries


Previous1Next

☒ one sample per file


Table 2. Calibration coefficients

coefficient	value
intercept	9.618237
coef1	-0.542061
coef2	0.000000
coef3	0.000000

Zoom OutHelp

Plot 1. GPC chromatogram of standards

x = 6.2717, y = 0.9971Add RT

Plot 2. GPC calibration curve

The second tab "Data" is used to load GPC files and display the results (**Table 3 Results**).

Calibration:

linear

polynomial 3

Enter the Intercept

10.99

and the Slope

-0.65

for a polynomial calibration curve, enter the 2nd

0

and 3rd coefficients

Calibration

Data

Select GPC file(s)

Browse...

2 files

Help

Table 3. Results

	FileName	peak	PDI	Mn	Mw	Mp	time1	time2	time_max	Af	Tf
1	data3.arw	1	1.15	5485	6315	6735	8.15	9.73	8.62	1.65	1.33
2	data5.arw	1	1.14	161311	183567	196804	6.28	7.68	6.68	1.68	1.36
3	data5.arw	2	1.08	33792	36481	40899	7.33	8.23	7.58	2.45	1.83

Showing 1 to 3 of 3 entries

Previous


1

Next

☐ Add calibration information

Save

It shows chromatograms and molecular weight distribution (MWD) graphs.



Calibration

There are two ways to set up calibration.

Calibration:

☒ linear
☐ polynomial 3

Enter the Intercept

10.99

and the Slope

-0.65

for a polynomial calibration curve, enter the 2nd

0

and 3rd coefficients

0

☒ Use these coefficients

1. If you know the *intercept* and *slope* coefficients, you can enter these values in the corresponding fields in the sidebar. Then check the **Use these coefficients** checkbox.

For linear calibration:

- set the Intercept (it is greater than 0);
- the Slope (it is negative number);
- leave the second and the third coefficients equal to 0.

For a polynomial fit:

- switch the **Calibration** to "**polynomial 3**";
- enter *intercept* and other coefficients;
- note that the coefficients must be very precise; **do not round them for polynomial fit.**

2. The second way is to load GPS data files with calibration standards on the "Calibration" tab.

The application does NOT read raw .dat files. You need to export .dat files to .txt, .arw, excel or .csv/.tsv files.

Typical GPC files have tab-separated format with some metadata at the beginning of the file. If there are multiple detectors, the application will only use data from the first one.

Select calibration file(s)

Browse... No file selected

sampleNarw - Notepad

```
File Edit Format View Help
"Sample Set Id" 1752
"Sample Set Name" "sad100124HPA"
"System Name" "System III"
"Instrument Method Name" "System 3_RI"
"SampleName" "HPAbHEA1"
"Injection" 1
"Vial" "12"
"Channel1" "410 "
0.0166667 -0.01404255
0.0333333 -0.004680851
0.05 0
0.0666667 0
0.0833333 0
0.1 0
0.116667 0
0.133333 0
0.15 0
0.166667 0
0.183333 0
```

sampleGarw - Notepad

```
File Edit Format View Help
"SampleName" "Vial" "System Name" "Date Acquired" "Injection"
"19NVPVPEMDia" "21" "System III" "8/3/2022 1:00:03 AM EDT"
0.0166667 0.01404255
0.0333333 0.004680851
0.05 0
0.0666667 0
0.0833333 0
0.1 0
0.116667 0
0.133333 0
0.15 0
0.166667 0
0.183333 0
0.2 0
0.216667 0.004680851
0.233333 0.004680851
0.25 0.009361702
0.266667 0.009361702
0.283333 0.009361702
0.3 0.009361702
```

PMMA12_8kDa.txt - Notepad

```
File Edit Format View Help

<CHANNEL_INFORMATION>
Channel ID Detector type Detector name Detector units Chan
1 RI RID: RI Signal nRIU
2 UV 260nm VwD: Signal A, 260 nm mAU
</CHANNEL_INFORMATION>

<RAW_DATA>
RT (mins) Response Trace 1 Response Trace 2
0 994.34 -16.6875123977661
0.008333333333333333 999.682982978829 -24.1852805527464
0.0166666666666667 1005.06528674655 -31.1419743399819
0.025 1010.40036722278 -33.8845981270779
0.0333333333333333 1015.69835997851 -32.1068812341849
0.0416666666666667 1020.96680879189 -28.0542818720725
0.05 1026.23978019587 -24.1734229988324
0.0583333333333333 1031.37792536259 -22.1674057058394
0.0666666666666667 1036.43910500869 -22.3559780900293
```


There are two options.

2.1. Load **multiple files** at once if each calibration sample is in a separate file (**each calibration standard has been injected into the GPC separately**).

Application extracts Mp (peak molecular weight) values for each calibration sample from the **SampleName** or **FileName** (if there is no sample name).

Examples of the sample name and Mp (in Da units):

PMMA12.8kDa	12800
PMMA1100	1100
PMMA31000	31000
PS2.55K	2550
PS200k	200000



These values are listed in the second column of **Table 1**.

If these values are incorrect or missing, double-click the appropriate cells in the **Table 1** and enter the new Mp(s) values.

The retention time values are given in the third column of **Table 1**. If the peak intensity is too low, the peak may not be detected, in which case the RT is missing.

Table 1. Calibration samples

	FileName	Mp	RT
1	calib_1.arw	1100	12.3
2	calib_2.arw	2200	11.7667
3	calib_3.arw	5090	11.15
4	calib_4.arw	9680	10.75
5	calib_5.arw	31000	9.9667
6	calib_6.arw	62200	0
7	calib_7.arw	102000	9.1333
8	calib_8.arw	200000	8.8833

Showing 1 to 8 of 8 entries Previous 1 Next

There are missing Mp values for the calibration samples, shown as 0 in Table 1. These values will be excluded from the calculations. To change Mp or retention time, double-click the corresponding cells in Table 1, enter the new values

If necessary, you can also adjust the peak retention time values (double click on table 1 cells).

Any missing values or values equal to zero will not be used to calculate calibration coefficients.


2.2. If a mixture of calibration standards has been injected into the GPC, load one or more files with chromatograms having multiple peaks and uncheck the "one sample per file" checkbox.

If only one calibration file is provided, the application will automatically assume that all calibration samples are in this file and will try to find all peaks and shoulder peaks.

Add the Mp values to the 2nd column of **Table 1** (double-click the corresponding cells). If there are shoulder peaks, these are also included in **Table 1**. So if a false peak has been detected and you do not want it to be used for the calibration, simply leave the value for Mp blank or 0, or double-click the corresponding RT cells and delete the incorrect value

The retention time values can also be added from **Plot 1**: click on this plot, the x and y values will appear in the text box below. To add this RT value to the **Table1**, click the "**Add RT**" button and add the Mp value (double-click the table cell).

To zoom in on this graph, drag the mouse from top left to bottom right across the graph area. To zoom out click on "**Zoom Out**".



The calibration chromatogram can be trimmed by adding a value for the system peak. In some cases, this may improve the calculation, for example, if the analyte peak is small.

In the sidebar, enter a value in the "**System Peak, Calibration**" field. The application will automatically recalculate the calibration. However, if Mp values were missed, you will need to re-enter these values in **Table 1**.

3. Choose *linear* or *polynomial3* fit for the calibration curve.

Calibration:

☐ linear

☒ polynomial 3

Enter the Intercept

10.99

and the Slope

-0.65

for a polynomial calibration curve,
enter the 2nd

0

and 3rd coefficients

0

☐ Use these coefficients

Table 2. Calibration coefficients


coefficient	value
intercept	9.618237
coef1	-0.542061
coef2	0.000000
coef3	0.000000




Table 2. Calibration coefficients

coefficient	value
intercept	130.934249
coef1	-37.877438
coef2	3.802035
coef3	-0.128143

Plot 2. GPC calibration curve



Plot 2. GPC calibration curve



Data tab

Load the GPC file(s) in the second tab.

Select GPC file(s)

Browse...

No file selected

Results table

Results are summarized in the **Table 3**. For each GPC file (**FileName**) and each peak (**peak**) there are:

- calculated M_p , average molecular weights (M_n and M_w), and PDI ;
- **time1** and **time2** are the start and end times of the peak, respectively; these values can be changed on the sidebar (see *Peak start and end*);
- **time_max** is the retention time at the peak's maximum;
- asymmetry factor (**Af**) and tailing factor (**Tf**) are calculated as given in [2];
- Check **Add calibration information** box to add it to the **Table 3**

Table 3. Results

	FileName	peak	PDI	Mn	Mw	Mp	time1	time2	time_max	Af	Tf	intercept	coef1	coef2	cc
1	data3.arw	1	1.15	5485	6315	6735	8.15	9.73	8.62	1.65	1.33	10.36087	-0.75813	0	
2	data5.arw	1	1.14	161311	183567	196804	6.28	7.68	6.68	1.68	1.36	10.36087	-0.75813	0	
3	data5.arw	2	1.08	33792	36481	40899	7.33	8.23	7.58	2.45	1.83	10.36087	-0.75813	0	

Showing 1 to 3 of 3 entries

Previous

1

Next

☒ Add calibration information

Save

In the sidebar you can switch between **linear** and **polynomial** calibration. The results will be updated automatically.

Calibration:

☐ linear

☒ polynomial 3

Enter the Intercept

10.99

and the Slope

-0.65

for a polynomial calibration curve, enter the 2nd

0

and 3rd coefficients

Calibration Data

Select GPC file(s)

Browse... 2 files

Help

Table 3. Results

	FileName	peak	PDI	Mn	Mw	Mp	time1	time2	time_max	Af	Tf	intercept	coef1	coef2	coef
1	data3.arw	1	1.14	5720	6507	6957	8.15	9.73	8.62	1.65	1.33	22.0388	-5.16019	0.54523	-0.02
2	data5.arw	1	1.17	150814	175990	187321	6.28	7.68	6.68	1.68	1.36	22.0388	-5.16019	0.54523	-0.02
3	data5.arw	2	1.07	31314	33522	37376	7.33	8.23	7.58	2.45	1.83	22.0388	-5.16019	0.54523	-0.02

Showing 1 to 3 of 3 entries

☒ Add calibration information

Save


Previous 1 Next

Click on **"Save"** button to save results as .csv file.

Plots


There will be two types of plots for each GPC file.

- The first is the chromatogram, where the graph name matches the GPC file name and the peaks are labeled 1, 2, 3, etc. Short vertical lines marking the beginning and end of each peak.
- And the second type of plots is the molecular weight distribution.



Base line correction



GPC chromatograms may have baseline drift, shift, or other problems [3]. Therefore, the application performs intensity correction so that the peaks start at 0 intensity.



Baseline anomalies are more often present at the beginning of the chromatogram, and the application does not look for peaks in this part of the graph. Therefore, if the analyte peak starts close to time 0, it will not be detected.

System peak

System peaks may be present in GPC chromatograms [4]. If the chromatogram is noisy or the peak of interest is too small compared to the system peak, the application may not detect it correctly.



Note:

- the system peak cannot be set to less than 20% of the maximum retention time;
- there are two fields for entering the system peak value for the calibration chromatogram and for the samples.


System peak, calibration: 0

System peak, samples: 10

Peak start and end


The start and the end of the peak are marked with short vertical lines on the chromatograms.

To shift these boundaries, enter values between -20 and 35 (percentage of peak width) into the fields in the sidebar and check the **"shift start time"** (or/and **"shift end time"**) box to apply the changes.



Adding tails

If peaks are merged you can add tails: check **"Add tails"** box on the sidebar.




The application will try to fit the sigmoid equation given in [5] and [6] and add the missing part of the chromatogram. If a large portion of the GPC curve is missing, the results may be inaccurate. Note that this app does not make any peak broadening corrections for molecular weight averages and does not add Gaussian curve.

Smoothing

If the chromatogram is noisy, the application will apply smoothing before finding the peaks. However, in some cases the noise may be detected as peaks.

No peaks

If no peaks are detected, only the raw data will be plotted.



This application was created using R version **4.4.2** and R packages:

DT **0.33** purrr **1.0.2** dplyr **1.1.4.9000** stringr **1.5.1** tidyr **1.3.1** shiny **1.9.1**
here **1.0.1**

- [1] S. Mori and H. G. Barth, *Size Exclusion Chromatography*, 1st ed. Springer Berlin, Heidelberg, 1999.
- [2] R. Freitag and L. G. Berruex, "Separation and Purification of Biochemicals," in *Encyclopedia of Physical Science and Technology*, 3rd ed., Elsevier, 2003, pp. 651–673.
- [3] S. Kromidas, "Quantification in LC and GC: a practical guide to good chromatographic data," 2009.
- [4] S. Levin and E. Grushka, "System peaks in liquid chromatography: their origin, formation, and importance," *Anal. Chem.*, vol. 58, no. 8, pp. 1602–1607, 1986.
- [5] K. Heusser, R. Heusser, J. Jordan, V. Urech, A. Diedrich, and J. Tank, "Baroreflex Curve Fitting Using a WYSIWYG Boltzmann Sigmoidal Equation.," *Front. Neurosci.*, vol. 15, p. 697582, 2021.
- [6] H. Motulsky and A. Christopoulos, *Fitting Models to Biological Data Using Linear and Nonlinear Regression: A Practical Guide to Curve Fitting*. in Mathematics. Oxford University Press, 2004.